



NATIONS ENVIRONMENT PROGRAMME FD

16 September 1983

Determination of total dissolved cadmium in sea-water by differential pulse anodic stripping voltammetry

DRAFT-NOT TO BE CITE PROJET-NE PAS DIVULC

Reference Methods For Marine Pollution Studies No. 18

Prepared in co-operation with





This document has been prepared by the Intergovernmental Oceanographic Note: Commission (IOC) of UNESCO and the United Nations Environment Programme (UNEP) under projects FP/ME/0503-75-07, ME/0503-81-01 and FP/0503-77-03. The assistance of Dr. Jean-Marie Martin in the preparation of the document is gratefully acknowledged.

The document has been reviewed by the IOC Group of Experts on Methods, Intercalibration (GEMSI) whose comments have been and Standards integrated into the text.

For bibliographic purposes this document may be cited as:

UNEP/IOC: Determination of total dissolved cadmium in sea-water by differential pulse anodic stripping voltametry. Reference Methods for Marine Pollution Studies No. 18, UNEP 1983.





UNITED NATIONS ENVIRONMENT PROGRAMME

16 September 1983

Determination of total dissolved cadmium in sea-water by differential⁻ pulse anodic stripping voltammetry

Reference Methods For Marine Pollution Studies No. 18

Prepared in co-operation with



IOC

AL.

PREFACE

The Regional Seas Programme was initiated by UNEP in 1974. Since then the Governing Council of UNEP has repeatedly endorsed a regional approach to the control of marine pollution and the management of marine and coastal resources and has requested the development of regional action plans. The Regional Seas Programme at present includes ten regions and has over 120 coastal States participating in it. (1), (2)

One of the basic components of the action plans sponsored by UNEP in the framework of the Regional Seas Programme is the assessment of the state of the marine environment and of its resources, and of the sources and trends of the pollution, and the impact of pollution on human health, marine ecosystems and amenities. In order to assist those participating in this activity and to ensure that the data obtained through this assessment can be compared on a world-wide basis and thus contribute to the Global Environment Monitoring System (GEMS) of UNEP, a set of Reference Methods and Guidelines for marine pollution studies are being developed and are recommended to be adopted by Governments participating in the Regional Seas Programme.

The methods and guidelines are prepared in co-operation with the relevant specialized bodies of the United Nations system as well as other organizations and are tested by a number of experts competent in the field relevant to the methods described.

In the description of the methods and guidelines the style used by the International Organization for Standardization (ISO) is followed as closely as possible.

The methods and guidelines, as published in UNEP's series of Reference Methods for Marine Pollution Studies, are not considered as final. They are planned to be periodically revised taking into account the development of our understanding of the problems, of analytical instrumentation and the actual need of the users. In order to facilitate these revisions the users are invited to convey their comments and suggestions to:

International Laboratory of Marine Radioactivity International Atomic Energy Agency c/o Musee Oceanographique MC98000 MONACO

which is responsible for the technical co-ordination of the development, testing and intercalibration of Reference Methods.

- (1) UNEP: Achievements and planned development of the UNEP's Regional Seas Programme and comparable programmes sponsored by other bodies. UNEP Regional Seas Reports and Studies No. 1 UNEP, 1982.
- (2) P. HULM: A Strategy for the Seas. The Regional Seas Programme: Past and Future UNEP, 1983.

1. 2.	saved on refm	nethod18 on refmethods: 16 September 19	83
3.			
4.			
5.	DE	TERMINATION OF TOTAL DISSOLVED CADMIUM	IN SEA-WATER
6.			
7.		BY DIFFERENTIAL PULSE ANODIC STRIPPING	VOLTAMETRY
8.			
9. 10.			
10.		CONTENTS	
12.			
13.			Page
14.			
15.			
16.	1.	Scope and field of application	1
17.			. 1
18.	2.	References	Ŧ
19.	7	Opingiala	1
20. 21.	3.	Principles	_
22.	4.	Reagents	2
23.			
24.	5.	Apparatus	2
25.			
26.	6.	Cleaning procedure	4
27.	_		5
28.	7.	Sampling and pretreatment	,
29. 30.	8.	Analysia	6
31.	ο.	Analysis	
32.	9.	Reporting of results	7
33.	- •		
34.	10.	Estimation of precision and accuracy	7
35.			
36.			

.

.

6

-

.

8.	1. SCOPE AND FIELD OF APPLICATION
9.	
0. 1. 2. 3. 4. 5.	This reference method describes an electrochemical technique for the determination of total dissolved cadmium in sea and estuarine waters. Detection limit is 1 ng of cadmium per kilogram of sea-water using a mercury film rotating electrode (MFE).
6.	
7.	
8. 9.	2. REFERENCES
9. 0.	- · · · · ·
1. 2. 3. 4. 5.	MART, L. (1979) Prevention of Contamination and other Accuracy Risks in Voltametric Trace Metal Analysis of Natural Waters. Part 1 : Preparation steps, filtration and storage of water samples. Frezenius Z. Anal. Chem. 296, 350-357.
6. 7. 8. 9. 0.	MART, L., NURNBERG, H. W. and VALENTA, P. (1980) Prevention of Contamination and other Accuracy Risks in Voltametric Trace Metal Analysis of Natural Waters. Part III : Voltametric ultratrace analysis with a multicell designed for clean bench working. Frezenius Z. Anal. Chem. 300 - 362.
2. 3. 4. 5. 6. 7.	SCHAULE and PATTERSON, C. C. (1980) Lead in the marine environment. Proceedings of the International Experts Discussion on Lead Occurrence, Fate and Pollution in the Marine Environment, Rovinj, Yugoslavia, 18 - 22 October 1977. M. Branica and Z. Konrad (eds.), Pergamon Press, Oxford. pp. 31-43.
• • •	EGG. PRINCETOWN APPLIED RESEARCH (1979) Model 303 static mercury drop electrode operating and service manual.
2.	
3.	3. PRINCIPLES
+ •	
•	After compliant and contractments the envelopic placed in the
5. 7.	After sampling and pretreatment, the sample is placed in the voltametric cell. Cadmium is concentrated in a thin mercury film by
•	application of a sufficiently negative potential at the MFE. The stripping
•	involves oxidation of cadmium into solution by application of an increasing
•	potential with superimposed pulse and the resulting current is recorded.
•	

- 1 -

.

`

.

	4. REAGENTS
	For the analysis, use only demineralized water and reagents of
recog	nized analytical quality with as low as possible Cd concentrations.
4.1 MILLI	Demineralized water supply from a MILLI-RO-MILLI Q device from PORE S.A. or tridistilled.
4.2	Nitric acid (distilled or Merck suprapur).
4.3	Nitric acid (analytical reagent grade).
4.4	Hydrochloric acid (distilled or Merck suprapur).
4.5	Hydrochloric acid (analytical reagent grade).
	Saturated potassium chloride solution prepared from Merck suprapur
potas	sium chloride.
4.7	Mercury (hexadistilled or Merck suprapur).
4.8 conta (4.2)	Mercuric nitrate solution. Prepare a Hg(NO ₃) ₂ solution ining 5 mg/l of Hg ⁺⁺ ; weigh 5 mg of mercury (4.7); add HNO ₃ until complete dissolution and dilute with distilled water to 1 1.
4.9	Cadmium stock solution. 1 g/l Merck titrisol in 0.1 M HCl (4.4).
stand	
	Cadmium standard solution. Prepare daily a series of working ard solutions by appropriate dilution of the cadmium stock solution with demineralized water adding 10 ml of HCl (4.4) per litre.
	ard solutions by appropriate dilution of the cadmium stock solution with demineralized water adding IO ml of HCl (4.4) per litre.
	ard solutions by appropriate dilution of the cadmium stock solution
the (ard solutions by appropriate dilution of the cadmium stock solution with demineralized water adding IO ml of HCl (4.4) per litre. NOTE: The concentration of the Cd standard solutions depends on
the (ard solutions by appropriate dilution of the cadmium stock solution with demineralized water adding 10 ml of HCl (4.4) per litre. NOTE: The concentration of the Cd standard solutions depends on d levels anticipated in the samples to be analysed.
the (ard solutions by appropriate dilution of the cadmium stock solution with demineralized water adding 10 ml of HCl (4.4) per litre. NOTE: The concentration of the Cd standard solutions depends on d levels anticipated in the samples to be analysed. Ultrapure nitrogen. Al ₂ O ₃ (0.3u) available from Methrom, OSI, etc.
the (ard solutions by appropriate dilution of the cadmium stock solution with demineralized water adding 10 ml of HCl (4.4) per litre. NOTE: The concentration of the Cd standard solutions depends on d levels anticipated in the samples to be analysed. Ultrapure nitrogen.
the (4.11	ard solutions by appropriate dilution of the cadmium stock solution with demineralized water adding 10 ml of HCl (4.4) per litre. NOTE: The concentration of the Cd standard solutions depends on d levels anticipated in the samples to be analysed. Ultrapure nitrogen. Al ₂ O ₃ (0.3u) available from Methrom, OSI, etc.
the (4.11 4.12	ard solutions by appropriate dilution of the cadmium stock solution with demineralized water adding 10 ml of HCl (4.4) per litre. NOTE: The concentration of the Cd standard solutions depends on d levels anticipated in the samples to be analysed. Ultrapure nitrogen. Al ₂ O ₃ (0.3u) available from Methrom, OSI, etc. 5. APPARATUS <u>Sampling equipment</u> 5.1.1 Sampling bottles without rubber or metallic parts. Modific
the (4.11 4.12	ard solutions by appropriate dilution of the cadmium stock solution with demineralized water adding IO ml of HCl (4.4) per litre. NOTE: The concentration of the Cd standard solutions depends on d levels anticipated in the samples to be analysed. Ultrapure nitrogen. Al ₂ O ₃ (0.3u) available from Methrom, OSI, etc. 5. APPARATUS <u>Sampling equipment</u>


5.1.3 Polyethylene bags. 135. 136. 137. 5.2 Filtration equipment 138. 5.2.1 A laminar air flow bench. 139. 140. A filtration set. The polycarbonate 600 ml MILLIPORE under 141. 5.2.2 pressure filtration set can be used. Glass and pyrex are prohibited. 142. 143. 0.45 u membrane filters (Millipore or Nuclepore). 144. 5.2.3 145. 5.2.4 Micropipettes. 146. 147. 5.3 148. Analysis equipment 149. A vertical laminar air flow bench. 150. 5.3.1 151. A balance (200 g with a precision of + 0.001 g). 152. 5.3.2 153. 154. 5.3.3 Plastic gloves. 155. 156. 5.3.4 Micropipettes. 157. 5.3.5 Quartz and teflon beakers. 158. 159. 160. 5.3.6 An electrochemical analyser equipped with differential pulse mode (e.g. the P.A.R.) 174 A from PRINCETOWN APPLIED RESEARCH, METHROM 161. 162. analyser). 163. 164. A 50 ml teflon polargraphic cell. 5.3.7 165. 166. 5.3.8 A working electrode. A rotating mercury film on glassy carbon electrode (MFE) (e.g. the METHROM E 628 has proved to give 167. 168. acceptable results). This electrode can be constructed in the laboratory 169. (MART et al. 1980). 170. 171. 5.3.9 An Ag/AgC1 reference electrode (e.g. METHROM EA 427). 172. 173. 5.3.10 A Pt counter electrode (e.g. METHROM EA 202). 174. 175. 5.3.11 An X-Y recorder. 176. 177. 5.3.12 Ultrapure nitrogen. Nitrogen N 60 from l'AIR LIQUIDE 178. containing less than 0.1 ppm of 0_2 can be used without column for 179. removing 0, traces. In this case only a scrubbing bottle containing 180. water acidified to pHl with HCl (4.4) is used. If another type of nitrogen 181. is used an oxygen scrubbing system is necessary. It is constituted of Vanadium II solution and zinc amalgam for Vanadium II regeneration. For 182. 183. the preparation of this solution see EGG PRINCETOWN APPLIED RESEARCH 184. (1979).

186. 5.3.13 A chronometer. 187. 188. 5.3.14 A high pressure UV lamp (e.g. HANAU 150 W). 189. 190. 191. 192. 6. CLEANING PROCEDURE 193. 194. 195. All quartzware, bottles, teflon, beakers and cells must be cleaned 196. before being used. 197. 198. Use only demineralized water. All manipulations have to be carried 199. out by hands covered with plastic gloves. 200. 201. 6.1 Polyethylene bottles and filtration sets 202.1 203. Fill the bottles with 10% HCl (4.5) and heat to 70 C for 2 days. 204. Repeat this operation twice. The procedure is then repeated with 2% HCl (4.4). The bottles are transferred to the laminar air flow bench and 205. filled with 1% HCl (4.4) and heated to 70 C for 4 days. After that time, 206. the acid is decanted and the bottle filled with demineralized water 207. 208. acidified by adding 1 ml HCl (4.4) to 1 l of water. Each capped bottle is 209. placed in a polyethylene bag (5.13). 210. 211. 6.2 Teflon beakers and quartz items 212. 213. The same procedure as 6.1 can be used, but replace HCl with HNO3. 214. 215. 6.3 Filters 216. 217. The filters are soaked in cold 50% HCl (4.5) for 2 days. The acid is changed and the filters are left in a fresh solution of the same acid for 2 218. days. After this, they are rinsed several times with demineralized water 219. and soaked in 10% HCl (4.4) for a week. After a repeated rinsing with 220. ultrapure water, the filters are kept in 1% HCl (4.4) for longer storage. 221. 222. Several days before use the filters are rinsed with demineralized water and 223. conditioned (i.e. soaked with sea-water). 224. 225. 6.4 Pipette tips 226. 227. Several hundred pipette tips are heated to 80 C in 20% HCl (4.5) for 228. two days in a 2 1 Erlenmeyer flask. This procedure is repeated twice, each 229. time followed by a careful rinsing with ultrapure water. A final cleaning 230. and heating is performed in 10% HCl (4.4). After rinsing several times. 231. the pipette tips are dried and sealed into cleaned polyethylene bags (20 per 232.

bag).

- 4 -

234.	6.5 Polyethylene bags				
235.					
236. 237.	Polyethylene bags are first rinsed with normal tap water, then they are filled with 10% HCl (4.5) and immersed into the same acid in a				
238.	polyethylene container for one half day. After rinsing with dermineralized				
239.	water, the procedure is repeated, with 10% HC1 (4.4) for 1 day. The bags				
248.	are rinsed with demineralized water and dried in the laminar air flow				
241.	bench.				
242.					
243.					
244.					
245.	7. SAMPLING AND PRETREATMENT				
246.	71 SAMEING AND INCINCAMENT				
247.	<u> </u>				
247.	7.1 Sampling				
243.	7.1 <u>Samping</u>				
247. 250.	If possible use a small boat. The best technique is to move the boat				
251.	upstream and sample in front of the boat. The sample is transferred to a				
252.	500 ml polyethylene bottle (5.12) which in turn is placed in a clean				
253.	polyethylene bag (5.13).				
255.	polyechylene bag (5.15).				
255.					
255. 256.	7.2 Pretreatment				
257.	The samples must be filtered as soon as possible on 0.45 um filters				
258.	in a laminar air flow bench in the laboratory. After filtration, the				
259.	samples are acidified to 1% with HCl (4.4).				
260.					
261.	NOTE: Problems could arise with water samples from polluted coastal				
262.	areas with high organic content. In this case, the samples must be UV				
263.	irradiated prior to analysis. The procedure is described by MART et al.				
264.	1980. The UV irradiation device is represented in figure 1.				



Reflector

U.V. Lamp

Voltametric cell filled with sample

Quartz beaker

Aluminium foil

Distance block from teflon

Glass dish with water for separating the cell from the outer atmosphere

Depending on its origin and its dissolved organic matter content the 269. 270. sample is irradiated for 1 to 12 hours. 271. 272. 273. 274. 8. ANALYSIS 275. 276. 277. Weigh 50 g of the sample directly in the voltametric cell on the balance placed in the laminar air flow bench. Add 20 ul of HCl (4.4); 50 278. ul of a solution of HgNO₃ (4.8). The sample is purged with ultrapure 279. 280. nitrogen (4.10) for at léast 10 minutes. 281. 282. Working conditions: Electrolysis time 6, 10, 12 minutes. 283. Electrolysis time depends on Cd concentration of the sample to be analysed. 284. For unknown samples an exploratory run with a deposition time of 5 minutes 285. has to be performed. From this run the appropriate deposition time from 3 286. to 12 minutes can be deduced. 287. 288. Initial potential = -1 V vs Af/AgC1. This rather negative potential 289. is chosen to obtain homogeneous mercury film. 290. 291. At least two standard additions must be made because of possible 292. matrix effects. 293. 294. Between each measurement the MFE has to be polished with 295. Al₂0₂ (0.3u). Aluminium powder (4.12) is spread over a wet filter 296. paper and the polish is obtained by rotating the electrode in contact with 297. the powder. After polishing, rinse the MFE with demineralized water. A 298. final rinsing step is performed by rotating the electrode for several 299. minutes in demineralized water acidified to pH 1 with HNO, (4.2), 300. followed by a second rinse with demineralized water. 301. 302. Then a blank must be run. For this purpose demineralized water spiked with 20 ul HCl (4.4), 50 ul of HgNO, (4.8), 100 ul of KCl (4.6) 303. is placed in the voltametric cell. The deposition time is 6 minutes. The 384. 305. voltamogram is registered. If the blank value is below 1 to 5 ng/kg the 306. demineralized water can be replaced by the sample after having allowed the 307. electrode to rotate during 1 minute at 0.5 V vs Af/AgC1 in order to remove 308. the film. Between the standard addition let also the electrode rotate at 309. the same potential. 310. 311. Another alternative to this procedure is proposed by MART et al. (1980): 312. if working conditions are good, no background signal is measured and there 313. is no need to compensate it. 314. 315.

9. REPORTING OF RESULTS 317. 318. 319. $h_1 =$ height of the peak corresponding to the sample. 320. 321. h₂ = height of the peak corresponding to the sample spiked 322. with A ng of Cd. 323. 324. h_3 = height of the peak corresponding to the sample spiked 325. with 2A ng of Cd. 326. 327. 328. $h_2 - h_1$ and $h_3 - h_2$ should be equal. If there is a slight 329. difference, take: 330. 331. = $\frac{1+2}{2}$ to calculate the cadmium content (Xng) of the sample. 332. 333. 334. 335. The Cd content is: 336. 337. $X_{ng} = A_{ng} \times h_1$ 338. 339. The Cd concentration (ng/kg) of the sample is: 340. 341. $C = \frac{Xng}{M}$ 342. 343. 344. Where M is the weight (kg) of sample placed in the voltametric cell. 345. 346. NOTE: This calculation is only for linear calibration curves. 347. 348. 349. 350. 10. ESTIMATION OF PRECISION AND ACCURACY 351. 352. 353. 354. 10.1 Precision 355. Estimate precision (C.V. = coefficient of variation) of the 356. analytical procedure by taking 7 different subsamples from the original 357. sample. The precision must be estimated for different ranges of 358. concentrations, e.g. 1 ng/kg to 10 ng/kg, 10 ng/kg to 100 ng/kg, 100 ng/kg 359. to 1000 ng/kg. 360. 361. 362. 10.2 Accuracy 363. Participate to intercalibration exercises and/or analyse a known 364. certified standard of a matrix similar to the material under study. 365. 366. 367.

- 7 -

LIST OF REFERENCE METHODS FOR MARINE POLLUTION STUDIES

LISTE DES METHODES DE REFERENCE POUR LES ETUDES DE POLLUTION MARINE

- UNEP/WHD : Guidelines for monitoring the quality of coastal recreational and shellfish-growing waters. (Draft) Reference Methods for Marine Pollution Studies No. 1, UNEP 1984.
- UNEP/WHO : Determination of total coliforms in sea-water by the membrane filtration culture method. Reference Methods for Marine Pollution Studies No. 2 Rev. 1, UNEP 1983.
- PNUE/QMS : Détermination des coliformes totaux dans l'eau de mer par la méthode de culture sur membranes filtrantes. Méthodes de Références pour les Études de Pollution Marine No 2, Rév. 1. PNUE 1983.
- UNEP/WHO : Determination of faecal coliforms in sea-water by the membrane filtration culture method. Reference Methods for Marine Pollution Studies No. 3 Rev. 1, UNEP 1983.
- PNUE/OMS : Détermination des coliformes fécaux dans l'eau de mer par la méthode de culture sur membranes filtrantes. Méthodes de Références pour les Etudes de Pollution Marine No 3, Rév. 1, PNUE 1983.
- UNEP/WHO : Determination of faecal streptococci in sea-water by the membrane filtration culture method. Reference Methods for Marine Pollution Studies No. 4 Rev. 1, UNEP 1983.
- PNUE/OMS : Détermination des streptocoques fécaux dans l'eau de mer par la méthode de culture sur membranes filtrantes. Méthodes de Références pour les Etudes de Pollution Marine No 4, Rév. 1, PNUE 1983.
- UNEP/WHO : Determination of faecal colliforms in bivalves by multiple test tube method. Reference Methods for Marine Pollution Studies No. 5 Rev. 1, UNEP 1983.
- PNUE/OMS : Détermination des coliformes fécaux dans les bivalves par le test des tubes multiples. Méthodes de Références pour les Etudes de Pollution Marine No 5, Rév. 1, PNUE 1983.
- UNEP/FAO/IAEA : Guidelines for monitoring chemical contaminants in marine organisms. Reference Methods for Marine Pollution Studies No. 6, UNEP. (in preparation)
- UNEP/FAG/IOC/ Sampling of selected marine organisms and sample preparation for trace metal analysis. IAEA : Reference Methods for Marine Pollution Studies No. 7 Rev. 2, UNEP 1984.
- UNEP/FAO/IOC/ Determination of total mercury in selected marine organisms by cold vapour atomic IAEA : absorption spectrophotometry. Reference Methods for Marine Pollution Studies No. 3 Rev. 1, UNEP 1984.
- UNEP/FAO/IAEA : Determination of total arsenic in selected marine organisms by hydride generation atomic absorption spectrophotometry. (Draft) Reference Methods for Marine Pollution Studies No. 9, UNEP 1985.
- UNEP/FAO/IAEA : Determination of total selenium in selected marine organisms by hydride generation atomic absorption spectrophotometry. Reference Methods for Marine Pollution Studies No. 10, UNEP 1984.
- UNEP/FAO/IDC/ Determination of total cadmium, zinc, lead and copper in selected marine organisms by IAEA : flameless atomic absorption spectrophotometry. Reference Methods for Marine Pollution Studies No. 11 Rev. 1, UNEP 1984.
- UNEP/FAO/IAEA : Sampling of selected marine organisms and sample preparation for the analysis of chlorinated hydrocarbons. Reference Methods for Marine Pollution Studies No. 12 Rev. 1, UNEP 1984.
- UNEP/FAO/IAEA : Determination of methylmercury in selected marine organisms by gas chromatography. Reference Methods for Marine Pollution Studies No. 13, UNEP 1984.
- UNEP/FAO/IOC/ Determination of CDTs and PCBs in selected marine organisms by packed column gas IAEA : chromatography. Reference Methods for Marine Pollution Studies No. 14 Rev. 1, UNEP 1985.

- UNEP/IDC/IAEA : Monitoring of tar on marine beaches. (Draft) Reference Methods for Marine Pollution Studies No. 15, UNEP 1985.
- UNEP/IAEA : Determination of DDIs, PCBs, PCCs and other hydrocarbons in sea-water by gas chrometography. (Draft) Reference Methods for Marine Pollution Studies No. 16, UNEP 1982.
- UNEP/IAEA : Determination of DDTs, PCBs and other hydrocarbons in marine sediments by gas-liquid chromatography. (Draft) Reference Methods for Marine Pollution Studies No. 17, UNEP 1982.
- UNEP/IOC : Determination of total dissolved cadmium in sea-water by differential pulse anodic stripping voltammetry. (Draft) Reference Methods for Marine Pollution Studies No. 18, UNEP 1983.
- UNEP/IOC/IAEA : Determination of mercury in estuarine waters and suspended sediment by cold vapour atomic absorption spectrophotometry. (Draft) Reference Methods for Marine Pollution Studies No. 19, UNEP 1985.
- UNEP/IOC/IAEA : Monitoring of petroleum hydrocarbons in sediments. Reference Methods for Marine Pollution Studies No. 20, UNEP. (in preparation)
- UNEP/WHO/IAEA : Determination of total coliforms in sea-water by multiple test tube (MPN) method. (Draft) Reference Methods for Marine Pollution Studies No. 21, UNEP 1985.
- UNEP/WHO/IAEA : Determination of faecal coliforms in sea-water by multiple test tube (MPN) method. (Draft) Reference Methods for Marine Pollution Studies No. 22, UNEP 1985.
- UNEP/WHO/IAEA : Determination of faecal streptococci in sea-water by multiple test tube (MPN) method. (Draft) Reference Methods for Marine Pollution Studies No. 23, UNEP 1985.
- UNEP/WMO/IAEA : Sampling of aerosols and wet precipitation for analysis of chemical pollutants. (Draft) Reference Methods for Marine Pollution Studies No. 24, UNEP 1985.
- SPC/UNEP : Coral reef monitoring handbook. Reference Methods for Marine Pollution Studies No. 25, UNEP 1984.
- UNEP/IAEA : Determination of total mercury in marine sediments and suspended solids by cold vapour atomic absorption spectrophotometry. (Draft) Reference Methods for Marine Pollution Studies No. 26, UNEP 1985.
- UNEP/IAEA : Determination of total cadmium in marine sediments by flameless atomic absorption spectrophotometry. (Draft) Reference Methods for Marine Pollution Studies No. 27, UNEP 1985.
- UNEP/WHO/IAEA : Determination of staphylococcus aureus in sea-water and sewage by the membrane filtration culture method. Reference Methods for Marine Pollution Studies No. 28, UNEP. (in preparation)
- UNEP/WHO/IAEA : Determination of pseudomones aeruginosa in sea-water and sewage by the membrane filtration culture method. Reference Methods for Marine Pollution Studies No. 29, UNEP. (in preparation)
- UNEP/WHO/IAEA : Isolation/Enumeration of salmonella from sea-water and sewage. Reference Methods for Marine Pollution Studies No. 30, UNEP. (in preparation)
- UNEP/IAEA : Determination of total chromium in marine sediments by flameless atomic absorption spectrophotometry. (Draft) Reference Methods for Marine Pollution Studies No. 31, UNEP 1985.
- UNEP/IAEA : Determination of total cobalt in marine sediments by flameless atomic absorption spectrophotometry. (Draft) Reference Methods for Marine Pollution Studies No. 32, UNEP 1985.
- UNEP/IAEA : Determination of total copper in marine sediments by flameless atomic absorption spectrophotometry. (Draft) Reference Methods for Marine Pollution Studies No. 33, UNEP 1985.

UNEP/IAEA	:	Determination of total lead in marine sediments by flameless atomic absorption spectrophotometry. (Draft) Reference Methods for Marine Pollution Studies No. 34, UNEP 1985.
UNEP/IAEA	:	Determination of total nickel in marine sediments by flameless atomic absorption spectrophotometry. (Draft) Reference Methods for Marine Pollution Studies No. 35, UNEP 1985.
UNEP/IAEA	:	Determination of total vanadium in marine sediments by flameless atomic absorption spectrophotometry. (Draft) Reference Methods for Marine Pollution Studies No. 36, UNEP 1985.
UNEP	;	Sampling and identification of common Mediterranean Scyphomedusae and evaluation of their occurrence. (in preparation)
UNEP/ICC/IAEA	:	Monitoring of petroleum hydrocarbons in sea-water. (in preparation)
UNEP/IAEA	:	Guidelines for monitoring of estuarine waters and suspended matter. (in preparation)
UNEP/WHO/IAEA	:	Determination of faecal coliforms in estuarine waters, suspended matter and sediments. (in preparation)
UNEP/WHO/IAEA	:	Determination of phosphorus in suspended matter and sediments. (in preparation)
UNEP/WHO/IAEA	;	Determination of nitrogen in suspended matter and sediments. (in preparation)
UNEP/WHO/IAEA	:	Determination of BOD ₅ and COD in estuarine waters. (in preparation)
UNEP/FAO/IAEA	:	Acute toxicity tests. (in preparation)
UNEP/IOC/IAEA	:	Determination of total cadmium in estuarine waters and suspended matter. (in preparation)
UNEP/FAO/IAEA	:	Biological non-acute toxicity tests. (in preparation)
UNEP/IOC/IAEA	:	Determination of basic oceanographic and meteorological conditions. (in preparation)
UNEP/IOC/IAEA	:	Determination of standard physical and chemical parameters. (in preparation)
UNEP/WHO/IAEA	:	Statistical methods for the evaluation of results from monitoring the quality of coastal recreational and shellfish-growing waters. (in preparation)
UNEP/FAO/IOC/ IAEA :		Determination of DDTs and PCBs in selected marine organisms by capillary column gas chromatography. (in preparation)
UNEP/IAEA	:	Determination of selected trace metals in aerosol and in wet precipitation. (in preparation)
UNEP/IAEA	:	Determination of halogenated hydrocarbons in aerosol and in wet precipitation. (in preparation)
UNEP/WMO/IAEA	:	Sampling of dry deposition. (in preparation)
UNEP/WHO/IAEA	:	Determination of methylmercury, total mercury and selenium in human hair. (in preparation)
UNEP/WHO/IAEA	:	Guidelines for monitoring and epidemiological studies on health effects of methylmercury. (in preparation)
UNEP/IOC/IAEA	:	Guidelines for the determination of riverine inputs of contaminants to estuaries. (in oreparation)

Issued and printed by:

Programme Activity Centre for Oceans and Coastal Areas United Nations Environment Programme

Additional copies of this and other publications issued by the Programme Activity Centre for Oceans and Coastal Areas can be obtained from:

Programme Activity Centre for Oceans and Coastal Areas United Nations Environment Programme Palais des Nations GENEVA Switzerland

